

BLOOD SAMPLING (BLOOD COLLECTOR CARDS)

HOW TO COLLECT BLOOD SAMPLES USING BLOOD CARDS

Always wash your hands or use clean gloves. Allow blood samples to dry completely before shipping, but do not dry them in direct sunlight or in a plastic bag. When collecting blood, use a new needle, syringe or pin prick device for each animal. When collecting, packaging and submitting samples for DNA testing, it is critical to eliminate cross-contamination. Rinsing needles is not enough to prevent possible DNA residue. Please follow these simple steps to improve accuracy of test results:

1. Verify the animal's tattoo and/or tag number and clearly record it on the blood card in the space provided. This ID will be used for reporting results, so it's important to double-check it for accuracy.
2. Using a new needle or piercing device for each animal, prick a vein in the animal's ear. As blood forms at the site of the prick, open the blood card and touch it to the drop on the indicated square. Apply only enough blood to the card to fill the two squares outlined, but not so much that the card is soggy and overly saturated. Two or three heavy drops of blood per square is plenty. A syringe may also be used to collect and spot the blood cards.
3. Next, bend the top flap of the card up so that the sample dries while the card is in an open position.
4. Set the blood card aside with the top flap open until dry. Make sure that blood cards do not touch—this is a possible source of DNA cross-contamination. Also, do not dry blood cards in direct sunlight, as UV radiation destroys DNA.
5. Once dry, rubber band closed blood cards together in batches of 10 – 15 and combine several batches together into resealable plastic storage bags.
6. Record the Animal ID along with other animal information on the Sample Information Form, such that each sample can be correctly matched with each animal. Barcodes are assigned to the samples upon arrival in the lab. Samples may be shipped by regular or express mail to Zoetis.



All trademarks are the property of Zoetis Inc., its affiliates and/or its licensors. All other trademarks are the property of their respective owners. ©2013 Zoetis Inc. All rights reserved. ZAG13007

GLOBAL HEADQUARTERS

333 Portage Road, Building 300 • Kalamazoo, Michigan 49007-9970

Phone: 877-233-3362 • Fax: 269-833-4711

www.zoetis.com